TESLA SMART OXIMETER USER MANUAL

EN

DE

CZ

SK


PL

HU

SL

THIS MEDICAL INSTRUMENT
MUST BE USED ACCORDING TO
INSTRUCTIONS TO ENSURE
ACCURATE READINGS

1. Power button
2. Display
3. SpO₂ data
4. Bluetooth/Battery indicator
5. Pulse Rate data
6. Plethysmogram
7. Battery cover

1 SAFETY

1.1 Instructions for the Safe Operation and Use of the Tesla Smart Oximeter

- Do not attempt to service the Tesla Smart Oximeter. Only qualified service personnel should attempt any needed internal servicing.
- Prolonged use or the patient's condition may require changing the sensor placement periodically. Change the sensor placement and check sensor is in contact with exposed skin, circulatory status, and correct alignment at least every 2 hours.
- Oxygen saturation (SpO₂) measurements may be adversely affected in the presence of high ambient light. The sensor area should be shielded with a surgical towel, plaster, or bandage if necessary.
- The following factors may cause interference to the testing accuracy of the Tesla Smart Oximeter:

- High-frequency electrosurgical equipment.
- Placement of the sensor on an extremity with a blood pressure impacted by arterial catheter, or intravascular line.
- If the patient has severe hypotension, vasoconstriction, severe anemia, or hypothermia.
- If the patient is in cardiac arrest or in shock.
- Fingernail polish or false fingernails may cause inaccurate SpO₂ readings.
- The device should be given at least 10 minutes to come to normal working temperature if stored in a hot or cold environment.
- The device is non-sterile and not intended to be sterilized.

1.2 Warning

The MEDICAL ELECTRONIC EQUIPMENT is suitable for home healthcare environments:

- Although the medical electronic equipment conforms to the intent of the standard EN 60601-1-2 in relation to electromagnetic compatibility, electrical equipment may produce interference. If interference is suspected, move the equipment away from the sensitive device.
- Portable and mobile RF communication equipment can affect this instrument's normal operation.
- Explosion hazard - Do not use the Tesla Smart Oximeter in a flammable atmosphere where concentrations of flammable anesthetics or other materials may occur.
- Do not throw batteries in fire as this may cause them to explode.
- Do not attempt to recharge normal dry-cell batteries, they may leak and

may cause a fire or even explode.

- Do not use the Tesla Smart Oximeter in an MRI or CT environment.
- Do not modify this equipment without the authorization of the manufacturer.
- If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of equipment.
- Do not use near active high frequency surgical equipment and the RF shielded room of an medical electronic equipment system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.
- Use of this equipment adjacent to or stacked with other equipment should be avoided, because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers, and cables other than those specified or

provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Tesla Smart Oximeter, including cables specified by the manufacturer. Otherwise, it may result in a degradation of the performance of this equipment.

- High-pressure sterilization cannot be used on the device.

- IF ANY: a list of all cables and maximum lengths of cables (if applicable), transducers and other accessories that are replaceable by the responsible organization and that are likely to affect compliance with the requirements of clause 7 (emissions) and clause 8 (immunity). accessories may be specified

either generically (e.g., shielded cable, load impedance) or specifically (e.g., by manufacturer and equipment or type reference).

- **IF ANY:** the performance that was determined to be essential performance and a description of what the operator can expect if the essential performance is lost or degraded due to disturbances (the defined term “essential performance” need not be used).

Cautions:

- Keep the operating environment free of dust, vibrations, corrosive, or flammable materials, and extremes of temperature and humidity.
- Do not operate the unit, if it is damp or wet because of condensation or spills. Avoid using the equipment immediately after moving it from a cold environment to a warm, humid location.
- Never use sharp or pointed objects to operate the front-panel switches.

- The batteries must be taken out from the battery compartment if the device will not be used for a long time.
- The device shall only be used with the battery cover closed.
- The batteries must be disposed of properly, according to local regulation after their use.
- Keep the device away from children and pets to avoid swallowing.

2 INTRODUCTION

2.1 Intended use

The Tesla Smart Oximeter is a non-invasive device intended for spot checking of functional oxygen saturation of arterial hemoglobin (SpO₂) and pulse rate (PR). This portable device is suitable for non-professional adult patients of any color in clinical institutions and home environments.

2.2 Brief Device Description

The Tesla Smart Oximeter is based on digital technology. The device is intended for non-invasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO₂). Advanced DSP algorithm can minimize the influence of motion artifact and improve measurement accuracy of low perfusion.

The Tesla Smart Oximeter can be used to measure human SpO₂ and pulse rate through the finger. The product is suitable for private home-use, hospitals (including clinical use in internist/surgery, anesthesia, pediatrics etc.), social medical organizations, physical care in sports etc.

Note:

- DSP algorithm: Digital signal processor algorithm.
- Low Perfusion: In physiology, perfusion is the process of a body

delivering blood to a capillary bed in its biological tissue. Under the condition of low perfusion, the measurement of non-invasive saturation of pulse-blood oxygen is not accurate.

2.3 Product Features

- Lightweight for easy carrying and use.
 - Manually adjust the direction of interface.
 - Color OLED display, simultaneous display for testing value and plethysmography.
 - Visual & Sound reminder function.
- Real-time spot-checks.
- Wireless Bluetooth for data transmission.
 - Low Battery voltage indicator.
 - Automatic switch-off.
 - Includes two standard AAA 1.5V alkaline batteries.

Note:

- **Plethysmograph:** is an instrument for measuring changes in volume within an organ or whole body (usually resulting from fluctuations in the amount of blood or air it contains).

CAUTION: The device must not be used on children below 3 years of age, as it cannot be guaranteed to provide accurate results.

CAUTION: The Tesla Smart Oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.

CAUTION: A function tester cannot be used to assess the accuracy of the Tesla Smart Oximeter monitor or sensor. Clinical testing is used to establish the SpO₂ accuracy. The measured arterial SpO₂ value (SpO₂) of the sensor is compared to arterial hemoglobin oxygen (SaO₂) value, determined from

blood samples with a laboratory CO-oximeter. The accuracy of the sensors in comparison to the CO-oximeter samples measured over the SpO₂ range of 70-100%. Accuracy data is calculated using the root-mean-square (Arms value) for all subjects. Only about two-thirds of pulse oximeter measurements can be expected to fall within \pm Arms of the value measured by a CO-oximeter. A pulse simulator shall be used to assess pulse rate accuracy. The measured pulse rate is compared to the preset pulse rate value in the simulator. Accuracy data is calculated using the root-mean-square (Arms value) for all subjects.

2.4 Expected Service Life

Is dependent on:

- The expected service life of the medical electronic equipment;
- The expected service life of parts or accessories shipped with the medical electronic equipment;
- Where the shelf life is less than the expected service life, the shelf life of parts or accessories shipped with the medical electronic equipment;

3 INSTALLATION, SETUP AND OPERATION

3.1 OLED display parameter setting

When the device is in the measuring interface, press the direction button for 1 second in order to enter the menu page (figure 3.1.1 and figure 3.1.2). There are two sub-menus.

3.1.1 Reminder Setup

Press the direction button for 1 second and enter the Reminder Setup. User can adjust the setting through moving the “*” symbol to the Sound Reminder, Beep, Restore or Brightness.

- Sound reminder

Press the direction button for 1 second, move the “*” symbol to the “Sound Reminder”, long press the direction button to turn the alarm sound on/off. If the measured value exceeds the maximum or minimum value of SPO2 or PR, there will be an alarm sound, if this function is turned on.

- Beep

Press the direction button for 1 second, move the “*” symbol to the back of “Beep”, long press the direction button to turn the beep sound on/off. When the beep is turned on, the sound emitted during the test indicates the pulse rate.

- Bluetooth

Press the direction button for 1 second, move the “*” symbol to the back of Bluetooth, long press the direction button to turn it on/off.

- Restore

With the “*” symbol next to “Restore”, long press the direction button. This will change to “OK”, which restores the device to factory settings.

- Brightness

With the “*” symbol next to “Brightness”, long press the direction button to change the Brightness value on a scale from 1 to 5.

3.1.2 Limit Value Setting

With the “*” symbol next to “Reminder Setup”, long press the direction button until the “Limit Setup” menu appears. You can then press the direction button to select the items. Press the direction button for 1 second to change to the data you need.

On the “Limit Setup” menu page, with the “*” symbol next to “+/-”, press the direction button for 1 second to change the “+” to “-” or change the “-” to “+”.

When “+” shows on the right side, by pressing the direction button for 1 second and moving the “*” next to the Spo2 Hi or PR Hi setting, you can increase the value to a higher value (until it reaches the maximum).

When “-” shows on the right side, by press the direction button for 1 second and moving the “*” next to the Spo2 Lo or PR Lo value setting, you can reduce the value to a lower value (until it reaches the minimum).

3.2 Operation

3.2.1 Installing batteries


Place two AAA batteries into the battery compartment taking care to align the polarities, then close the cover.

WARNING: Do not attempt to recharge normal alkaline batteries, as they may leak and may cause a fire or even explode.

3.2.2 Turning on and applying the Tesla Smart Oximeter

Put one of your fingers into the rubber opening of the Tesla Smart Oximeter with nail facing upward, then release the clamp. Press the power button to turn the Tesla Smart Oximeter on. The oximeter will automatically turn off, if there is no finger in the device for more than 16 ± 2 seconds.

3.2.3 Read data from display screen

Display	Description	Show	Note
OLED Screen	The screen display can scroll through four directions with six different display modes by pressing the power button.	 <p>The 'Show' column displays six different screen layouts. Each layout shows SpO2 (98%), PR (65 bpm), and P1% (5.5). The layouts vary in the placement of these values and the inclusion of a bar graph or a waveform.</p>	<ol style="list-style-type: none"> 1. When the battery power is at its lowest level, the battery symbol will be shown, reminding users to replace the batteries. 2. The plethymogram can be considered accurate if the wave symbol is fluctuating regularly.

4 APP OPERATION

4.1 Preparation for use

- Tesla Smart is available for both iOS and Android OS. Search the name ,Tesla Smart' in App Store or Android Market or scan the QR-Code to download the App and register or login your account.
- Make sure that the smartphone is connected to the 2.4Ghz network and have Bluetooth & GPS location turned on.

4.2 Connect the device

- Open Tesla Smart app and click “+” sign to Add new device
- Select Health care - Tesla Smart Oximeter

- Please follow the instructions to add the Tesla Smart Oximeter step by step on your phone

4.3 Control device

- Control oximeter with APP

5 CLEANING AND DISINFECTION

5.1 Cleaning

Switch off the power and take out the batteries before cleaning. Keep the exterior surface of the device clean and free of dust and dirt. Clean the exterior surface (display screen included) of the unit with a soft dry cloth. Use 75% medical alcohol to clean the surface by applying a small amount with a dry cloth to avoid the alcohol entering the device.

5.2 Disinfection

Disinfect the machine after each use, if multiple patients use the machine, e.g., in a hospital.

Use 75% medical alcohol to clean the surface that was in contact with the patient.

CAUTION: Do not use strong solvent, e.g., acetone.

CAUTION: Never use an abrasive such as steel wool or metal polish.

CAUTION: Do not allow any liquid to enter the product and do not immerse any part of the device in liquid.

CAUTION: Avoid pouring liquid on the device while cleaning.

CAUTION: Do not leave cleaning solution on the surface of the device.

6 TROUBLESHOOTING AND MAINTENANCE

6.1 Maintenance

Replace the batteries in a timely manner, if the battery indication is

low. Clean the surface of the oximeter before it is used for diagnosis of patients. Remove the batteries from the battery compartment if the oximeter will not be operated for a long time. It is best to store the product where the ambient temperature is between -25°C and 55°C and humidity is 15% to 93%. Regular inspection is recommended in order to make sure that no obvious damage is present that may affect the safety and performance of the device. Do not expose the device to flammable substances, high or low temperatures or humidity levels outside those in the operation conditions.

Problem	Possible Reason	Resolution
Oxyhemoglobin or heart rate cannot be shown.	<ol style="list-style-type: none"> 1. Finger is not inserted correctly. 2. Patient's perfusion is too low to be measured. 	<ol style="list-style-type: none"> 1. Retry by inserting the finger thoroughly. 2. Try a few more times to make sure there is no problem with the product itself. Otherwise seek medical help for exact diagnosis.
Oxyhemoglobin or heart rate is unstable.	<ol style="list-style-type: none"> 1. Finger might not be inserted correctly. 2. Finger is trembling, or patient's body is moving. 	<ol style="list-style-type: none"> 1. Retry by inserting the finger thoroughly. 2. Try to help the patient keep calm and still.
Oxyhemoglobin or heart rate is outside of the standard range.	<ol style="list-style-type: none"> 1. Finger might not be inserted correctly. 2. Patient's SPO₂ & PR is abnormal. 	<ol style="list-style-type: none"> 1. Retry by inserting the finger thoroughly. 2. Seek medical help for further examination.
The oximeter cannot be turned on.	<ol style="list-style-type: none"> 1. Batteries may need replacement. 2. Batteries might be installed incorrectly. 3. The oximeter might be damaged. 	<ol style="list-style-type: none"> 1. Please replace batteries. 2. Please reinstall the batteries. 3. Please contact your point of purchase.
The screen suddenly turns off.	<ol style="list-style-type: none"> 1. The device automatically turns off, if there is no signal detected for more than 16 seconds. 2. Batteries may need replacement. 	<ol style="list-style-type: none"> 1. This is a normal process. 2. Replace the batteries.

7 SPECIFICATION

Device Name: Tesla Smart Oximeter

Dimension (L×W×H): (62 × 34 × 31) mm

Anti-electric Shock Type:

Internally powered equipment

Anti-electric Shock Equipment Degree: Type BF

EMC type: Group I Class B

Enclosure Degree of ingress protection: IP22

Internal Power: 2×AAA, 1,5V alkaline batteries

Power consumption: Below 45mA

Screen: 0.96" OLED

SpO2 Display: 35-100%

Pulse rate Display: 30-250 BPM

Resolution: SpO2: 1%

Heart rate: 1 BPM

Measure Accuracy: SpO2: ±3% (70%-100%);

Unspecified (<70%)

PR: ±2BPM

Data averaging and other signal processing: 8s

Data Update Period: 1s

Operating Environment:

Temperature: +5°C to +40°C

Humidity: 15% to 93%, non-condensing

Air Pressure: 70Kpa-106Kpa

Storage & Transport Environment:

Temperature: -25C to +55C

Humidity: 15% to 93%

non-condensing

Air Pressure: 70Kpa-106Kpa

8 CLINICAL SPO2 ACCURACY

The below table shows statistic distribution of an invasive controlled desaturation study, which guided by ISO80601-2-61, Annex EE, Guideline for evaluating and documenting SpO2 Accuracy in human subjects. The statistic distribution displayed the accuracy distribution between the range of 70% ~ 100%, which may be helpful to user.

Item	100-70%	100-90%	90-80%	80-70%
Bs	1.0	0.80	1.26	0.86
Sres	1.94	0.86	1.61	2.91
Arms	2.18	1.17	2	2.99

Electromagnetic emissions

Emissions test	Compliance
RF emissions -CISPR 11	Group 1
RF emissions -CISPR 11	Class B
Harmonic emissions -IEC 61000-3-2	Not application
Voltage fluctuations/ flicker emissions -IEC 61000-3-3	Not application

Electromagnetic Immunity		
Immunity Test	IEC 60601-1-2 -Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines ±2 kV input/output lines ±1 kV	Not applicable
Surge IEC 61000-4-5	line(s) to line(s)±1 kV.; line(s) to earth±2 kV. 100 kHz repetition frequency	Not applicable
Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle; At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% 1 cycle; And 70% 25/30 cycles Single phase: at 0; 0% 300 cycles	Not applicable
Power frequency magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	150KHz to 80MHz: 3Vrms, 6Vrms (in ISM and amateur radio bands); 80% Am at 1kHz	Not applicable
Radiated RF -IEC61000-4-3	10 V/m; 80 MHz 2,7 GHz; 80 % AM at 1 kHz	10 V/m; 80 MHz 2,7 GHz; 80 % AM at 1 kHz

Guidance and manufacturer's declaration Electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	Immunity test Level (V/m)
Radiated RF IEC 61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	385	380-390	TETRA 400	Pulse Modulation 18 Hz	1,8	0.3	27
	450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710 745 780	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	810 870 930	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	2	0.3	28
	1720 1845 1970	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3,4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
	2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240 5240 5785	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9

9 MANUFACTURER'S EMC DECLARATION

- 1) All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.
- 2) Guidance and manufacturer's declaration -Electromagnetic Emissions and Immunity.

INFORMATION ABOUT DISPOSAL AND RECYCLING

This product is marked with the symbol for separate collection. The product must be disposed of in accordance with the regulations for disposal of electrical and electronic equipment (Directive 2012/19/EU on waste electrical and electronic equipment). Disposal together with regular municipal waste is prohibited. Dispose of all electrical and electronic products in accordance with all local and European regulations at the designated collection points which hold the appropriate authorisation and certification in line with the local and legislative regulations. Correct disposal and recycling help to minimise impacts on the environment and human health. Further information regarding disposal can be obtained from the vendor, authorised service centre or local authorities.

EU DECLARATION OF CONFORMITY

Hereby, Tesla Global Limited declares that the radio equipment type A310B is in compliance with 2014/53/UE directives.

The full text of the EU declaration of conformity is available at the following internet address: tsl.sh/doc

Connectivity: Bluetooth

Frequency band: 2.412 - 2.472 MHz

Max. radio-frequency power (EIRP): < 20 dBm

0123

Shenzhen Aeon Technology Co., Ltd.

RM6H02, Block 27-29, Tianxia IC Industrial Park, Majialong, No.133 of Yiyuan road, Nantou Street, Nanshan District, 518052, Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 80, 20537 Hamburg, Germany

Manufactured for:
Tesla Global Limited
www.teslasmart.com

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